

Instalacija LIS barkod štampača

[Instalacija na Linux OS](#)

[Instalacija na Windows OS](#)

[Problemi prilikom instalacije na Windows OS](#)

[Nedostaje MSVCP110.dll](#)

[Štampač ne štampa nalepnice](#)

Instalacija na Linux OS


Prilikom instalacije štampača u Linux operativnom sistemu potrebno je:

- Da ste ulogovani kao administrator
- U web pretraživač upisati: localhost:631/admin


Zatim pratiti sledeće korake

The screenshot shows the CUPS administration interface. The 'Server' section is expanded, showing 'Server Settings' with a 'Change Settings' button. Red arrows and numbers 1, 2, and 3 indicate the steps: 1 points to the 'Advanced' link, 2 points to the 'Change Settings' button, and 3 points to the 'Add Printer' button.

Klikom na dugme “Change settings” biće potrebno da se unese korisničko ime i lozinka. Nakon što se klikne na “Add Printer” prelazi se na sledeći ekran.



Izabere se lokalni barkode štampač i klikne na dugme “Continue”.




Ime štampača mora da bude isti kao naziv koji se nalazi u našoj bazi.



Iz sekcije Make izabere se Raw i klikne na dugme “Add Printer”. I do kraja se sve potvrdi. Potrebno je restartovati računar.


Instalacija na Windows OS

Prilikom instalacije štampača u Windows operativnom sistemu potrebno je da se pokrene LIS_LabelPrinterService-Setup-v1.0.exe



LIS Server je adresa servera, unosi se u formatu kao što je prikazano na slici:
`tcp://ip_adresa_servera:6663`

Printer Name je naziv koji ćete dobiti od LIS tehničke podrške. Nakon što su uneti svi podaci klikne se na dugme "Next >".






Nakon što se instalacija završi poželjno je restartovati računar, da bi se servis pokrenuo.

Problemi prilikom instalacije na Windows OS

Nedostaje MSVCP110.dll



Ukoliko se desi da se prilikom instalacija pojavi greška kao što je prikazana na gornjoj slici potrebno je instalirati *Visual C++ Redistributable for Visual Studio 2012 Update 4* koji se nalazi na linku <https://www.microsoft.com/en-us/download/details.aspx?id=30679> .

Za verziju Windows 10, Microsoft Visual C++ Redistributable možete preuzeti na linku

<https://www.itechtics.com/microsoft-visual-c-redistributable-versions-direct-download-links/>

Štampač ne štampa nalepnice

Pre svega treba proveriti da li je servis LISLabelPrinter pokrenut na Start -> Control Panel -> Administrative Tools -> Services, kao što je prikazano na sledećoj slici.

Windows Explorer window showing the contents of the Administrative Tools folder. The address bar shows the path: << All Control Panel Items >> Administrative Tools. The search bar contains the text "Search Administrative Tools".

The left sidebar shows the following navigation options:

- Favorites: Desktop, Downloads, Recent Places
- Libraries: Documents, Music, Pictures, Videos
- Computer: Local Disk (C:), VirtualBox-Folders, Network


The main pane displays a list of Administrative Tools items in a table format:

Name	Date modified	Type	Size
Component Services	7/14/2009 6:57 AM	Shortcut	2 KB
Computer Management	7/14/2009 6:54 AM	Shortcut	2 KB
Data Sources (ODBC)	7/14/2009 6:53 AM	Shortcut	2 KB
Event Viewer	7/14/2009 6:54 AM	Shortcut	2 KB
iSCSI Initiator	7/14/2009 6:54 AM	Shortcut	2 KB
Local Security Policy	9/7/2018 9:35 PM	Shortcut	2 KB
Performance Monitor	7/14/2009 6:53 AM	Shortcut	2 KB
Print Management	9/7/2018 9:35 PM	Shortcut	2 KB
Services	7/14/2009 6:54 AM	Shortcut	2 KB
System Configuration	7/14/2009 6:53 AM	Shortcut	2 KB
Task Scheduler	7/14/2009 6:54 AM	Shortcut	2 KB
Windows Firewall with Advanced Security	7/14/2009 6:54 AM	Shortcut	2 KB
Windows Memory Diagnostic	7/14/2009 6:53 AM	Shortcut	2 KB
Windows PowerShell Modules	7/14/2009 7:32 AM	Shortcut	3 KB


The bottom status bar shows details for the selected "Services" item:

- Services Shortcut
- Date modified: 7/14/2009 6:54 AM
- Size: 1.25 KB
- Date created: 7/14/2009 6:54 AM


Kada su podešavanja usklađena sa podešavanjima LIS-a ukoliko štampač ne štampa nalepnice treba proveriti da li je Firewall isključen na računaru, a ako nije potrebno ga ja isključiti.



Potrebno je kliknuti na "Turn Windows Firewall on or off" i isključiti Firewall kao što je prikazano na sledećoj slici.



Nakon toga, potrebno je vratiti se na prethodni ekran i kliknuti na "Advanced settings".



I tu je potrebno ugasiti Firewall klikom na "Windows Firewall Properties" i podesiti podešavanja kao na sledećoj slici.

